



PRACTIC 2

Manual Hospital Bed 2 Cranks

Detaylı ürün bilgisi için
QR Kodu okutun.
Scan QR Code for detailed
product information.



Genel Özellikler

- Krank mekanizması ile sırt ve bacak ayarı
- Kolay gizlenebilir ve kilitlenebilir korkuluklar
- Kolay çıkarılıp takılabilen PP baş ve ayak panelleri
- Fowler ve vasküler auto-contour pozisyonlar
- Dört bölmeli çelik yatma yüzeyi platformu
- Nano seramik kaplı elektrostatik toz boyalı metal aksam
- Cırcır ile ayak dayanağı ayarı
- Yüksekliği ayarlanabilir serum askısı
- 125 mm çapında 360° dönebilen çapraz kilitlenebilir tekerler
- 4 köşede koruyucu tamponlar
- 4 köşede serum askısı ve aksesuar yuvaları

General Features

- Back and leg adjustment with crank mechanism
- Easily concealable, lockable side rails
- Easy to remove and install PP head and foot panels
- Fowler and vascular auto-contour Positions
- Four-compartment steel lying surface platform
- Nano ceramic coated electrostatic powder painted metal parts
- Ratchet footrest adjustment
- Height adjustable IV pole
- 125 mm diameter 360° swivel diagonally lockable castors
- Protective bumpers at all four corners
- IV Pole and accessory slots at the four corners



Uzunluk
Length
213 cm



Genişlik
Width
107 cm



Yükseklik (Max)
Height (Max)
50 cm



Yatak Platformu
Bed Platform
90 x 190 cm



Sırt Hareketi
Back Movement
70°



Bacak Dayama
Leg Rest Angle
40°



Teker Çapı
Castor Diameter
125 mm



Güvenli Çalışma Yüğü
Safe Working Load
200 kg

