



Relax L

Foldable Attendant Sofa Bed

Detaylı ürün bilgisi için
QR Kodu okutun.

Scan QR Code for detailed
product information.



Genel Özellikler

- Koltuk, kullanımı kolay manuel bir sistemle genişletilerek koltuğu bireysel bir yatağa dönüştürebilir
- Dirsek desteği olan kolçaklar yatak pozisyonu aldığı anda baş kısmını destekler
- Özel mekanizma yardımıyla sırt kısmı 6 kademeli olarak omurga açısına göre ayarlanabilir
- Yüksek yoğunluklu sünger yatak platformu sayesinde ergonomik olarak kullanılabilir
- Köpükler lateks içermeyen, nefes alabilen, temizlenmesi kolay, aşınmaya dayanıklı, mantar önleyici, antibakteriyel ve yanmaz deri ile kaplanmıştır

General Features

- The seat can be expanded with an easy to use manual system, turning the seat into an individual bed.
- Armrests with elbow support support the head when in bed position
- With the help of a special mechanism, the back can be adjusted in 6 stages according to the angle of the spine
- Can be used ergonomically thanks to the high density foam mattress platform
- The foams are covered with latex-free, breathable, easy to clean, abrasion resistant, anti-fungal, antibacterial and fireproof leather.



Uzunluk
Length
70 cm



Genişlik
Width
193 cm



Yükseklik (Max)
Height (Max)
46 cm



Yatma Pozisyonu
Lying Position
70x193 cm